Introduction to ALKfusiondiscover

2022-07-18

Introduction

This is an introduction to the **ALKfusiondiscover** package which can be used to evaluate whether EML4-ALK is present in blood plasma.

This package was created in order to increase the sensitivity of EML4-ALK detection from commercially available NGS products such the AVENIO (Roche) pipeline.

Paired-end sequencing of cfDNA generated BAM files can be used as input to discover EML4-ALK variants. This package was developed using position deduplicated BAM files generated with the AVENIO Oncology Analysis Software. These files are made using the AVENIO ctDNA surveillance kit and Illumina Nextseq 500 sequencing. This is a targeted hybridization NGS approach and includes ALK-specific but not EML4-specific probes.

The package includes six functions.

The output of first function, EML4_ALK_detection(), is used to determine whether EML4-ALK is detected and serves as input for the next four exploratory functions characterizing the EML4-ALK variant. The last function EML4_ALK_analysis() combines the output of the exploratory functions.

To serve as examples, this package includes BAM files representing the EML4-ALK positive cell line H3122 and the EML4-ALK negative cell line, HCC827.

Installation

Use devtools to install ALKfusiondiscover.

```
if (!require(devtools)) install.packages('devtools')
library(devtools)

install_github("CTrierMaansson/ALKfusiondiscover")
library(ALKfusiondiscover)
```

Package data

BAM files from the cell lines, H3122 and HCC827, are included in the package and can be used as examples to explore the functions.

```
H3122_bam <- system.file("extdata", "H3122_EML4.bam", package = "ALKfusiondiscover")

HCC827_bam <- system.file("extdata", "HCC827_EML4.bam", package = "ALKfusiondiscover")
```

EML4_ALK_detection()

This function looks for EML4-ALK mate pair reads in the BAM file.

Input:

file

The name of the file which the data are to be read from.

genome

character representing the reference genome. Can be either "hg38" or "hg19". Default = "hg38".

mates

interger, the minimum number EML4-ALK mate pairs needed to be detected in order to call a variant. Default = 2.

Output:

If EML4-ALK is detected a data.frame with soft-clipped reads representing EML4-ALK is returned. Otherwise "No EML4-ALK was detected" is returned.

Examples:

```
head(EML4_ALK_detection(file = H3122_bam, genome = "hg38", mates = 2))
                                   sequences
mate position cigar
#>
#> 1 29223691 42299657 94M2S
#> 2 29223375 42299657 94M2S
#> 3 29223479 42299657 94M2S
#> 4 29223686 42299657 94M2S
#> 5 29223636 42299657 94M2S
#> 6 29223687 42299657 94M2S
EML4_ALK_detection(file = HCC827_bam, genome = "hg38", mates = 2)
#> [1] "No EML4-ALK was detected"
```

EML4_sequence()

This function identifies the basepairs leading up to the EML4 breakpoint.

Input:

reads

data.frame returned by EML4 ALK detection().

basepairs

integer, number of basepairs identified from the EML4-ALK fusion. Default = 20.

Output:

If EML4-ALK is detected, returns a table of identified EML4 basepairs with the number of corresponding reads for each sequence. Otherwise "No EML4-ALK was detected" is returned.

Examples:

```
EML4_sequence(EML4_ALK_detection(file = H3122_bam, genome = "hg38", mates = 2), basepairs = 20)
#> EML4_seq
#> CCAGGCTGGAGTGCAGTGGT GGAGTGCAGTGGTGTATTT TCAGGCTGGAGTGCAGTGGT
#> 201 1 1
EML4_sequence(EML4_ALK_detection(file = HCC827_bam, genome = "hg38", mates = 2), basepairs = 20)
#> [1] "No EML4-ALK was detected"
```

ALK_sequence()

This function identifies the basepairs following the ALK breakpoint.

Input:

reads

data.frame returned by EML4 ALK detection().

basepairs

integer, number of basepairs identified from the EML4-ALK fusion. Default = 20.

Output:

If EML4-ALK is detected, returns a table of identified ALK basepairs with the number of corresponding reads for each sequence. Otherwise "No EML4-ALK was detected" is returned.

Examples:

```
ALK_sequence(EML4_ALK_detection(file = H3122_bam, genome = "hg38", mates = 2), basepairs = 20)
#> ALK_seq
#> CAGAATTTTAGCTTTGCAAT
                               CGGATTTTTAGCTTT CGGATTTTTAGCTTTTCATT
#>
                                             1
                                                                    2
                                            CT
#> CGGAATTTTAGCTTTGCATT
                                                                  CTG
#>
                                             8
                                                                    3
                       1
#>
                    CTGA
                                         CTGAT
                                                            CTGATTTT
#>
                      11
                                            16
                                                                    5
#>
               CTGATTTTT
                                    CTGATTTTTA
                                                         CTGATTTTTAG
#>
                                                                    3
                       6
#> CTGATTTTTAGATTTGCATT
                                  CTGATTTTTAGC
                                                       CTGATTTTTAGCT
#>
                                                                   10
                                            14
```

```
#>
         CTGATTTTTAGCTT
                              CTGATTTTTAGCTTT
                                                   CTGATTTTTAGCTTTG
#>
                                             3
                      10
#>
      CTGATTTTTAGCTTTGC
                           CTGATTTTTAGCTTTGCA
                                                CTGATTTTTAGCTTTGCAT
                       7
                                             8
#>
                                                                   1
#> CTGATTTTTAGCTTTGCATT CTGATTTTTAGCTTTGCAAT
                                                   CTGATTTTTAGCTTTT
#>
                      71
                                             1
                                                                   1
#> CTGATTTTTAGCTTTTCATA
                                  CTGATTTTTAT
                                                   CTGATTTTTATCTTTG
#>
                                             1
                                                                   2
#> CTGATTTTTATCTTTGCATT CTGATTTTTATCTTTTGATT CTGTGTTTTAGATTTGCATT
#>
                                             1
#> CTGTTTTTTTTTTTGCAAT
                                         CTGAA CTTATTTTTATCTTTGCATT
#>
                                             1
#>
              TTAGCTTTG
#>
ALK_sequence(EML4_ALK_detection(file = HCC827_bam, genome = "hg38", mates = 2), basepairs = 20)
#> [1] "No EML4-ALK was detected"
```

break_position()

This function identifies the genomic position in EML4 where the breakpoint has happened.

Input:

reads

data.frame returned by EML4_ALK_detection().

Output:

If EML4-ALK is detected, returns a table of genomic positions with the number of corresponding reads for each sequence. Otherwise "No EML4-ALK was detected" is returned.

Examples:

```
break_position(EML4_ALK_detection(file = H3122_bam, genome = "hg38", mates = 2))
#> break_pos
#> 42299750 42299757
#> 202    1
break_position(EML4_ALK_detection(file = HCC827_bam, genome = "hg38", mates = 2))
#> [1] "No EML4-ALK was detected"
```

break_position_depth()

This function identifies the read depth at the basepair before the breakpoint in EML4.

Input:

file

The name of the file which the data are to be read from.

reads

data.frame returned by EML4_ALK_detection().

Output:

If EML4-ALK is detected a single integer corresponding to the read depth at the breakpoint is returned. Otherwise "No EML4-ALK was detected" is returned

Examples:

```
break_position_depth(H3122_bam, EML4_ALK_detection(file = H3122_bam, genome = "hg38", mates = 2))
#> [1] 251
break_position_depth(HCC827_bam, EML4_ALK_detection(file = HCC827_bam, genome = "hg38", mates = 2))
#> [1] "No EML4-ALK was detected"
```

EML4_ALK_analysis()

This functions collects the results from the other functions of the package.

Input:

file

The name of the file which the data are to be read from.

genome

character representing the reference genome. Can be either "hg38" or "hg19". Default = "hg38".

mates

interger, the minimum number EML4-ALK mate pairs needed to be detected in order to call a variant. Default = 2.

basepairs

integer, number of basepairs identified from the EML4-ALK fusion. Default = 20.

Output:

A list object with clipped_reads corresponding to EML4_ALK_detection(), last_EML4 corresponding to EML4_sequence(), first_ALK corresponding to ALK_sequence(), breakpoint corresponding to break_position(), and read_depth corresponding to break_position_depth().

Examples:

```
H3122_results <- EML4_ALK_analysis(file = H3122_bam, genome = "hg38", mates = 2, basepairs = 20)
HCC827_results <- EML4_ALK_analysis(file = HCC827_bam, genome = "hg38", mates = 2, basepairs = 20)</pre>
head(H3122_results$clipped_reads)
                                                              sequences
mate position cigar
#> 1 29223691 42299657 94M2S
#> 2 29223375 42299657 94M2S
#> 3 29223479 42299657 94M2S
#> 4 29223686 42299657 94M2S
#> 5 29223636 42299657 94M2S
#> 6 29223687 42299657 94M2S
H3122_results$last_EML4
#> EML4_seq
#> CCAGGCTGGAGTGCAGTGGT GGAGTGCAGTGGTGTATTT TCAGGCTGGAGTGCAGTGGT
                             1
H3122_results$first_ALK
#> ALK_seq
#> CAGAATTTTAGCTTTGCAAT
                   CGGATTTTTAGCTTT CGGATTTTTTAGCTTTTCATT
              1
                             1
                                           2
#> CGGAATTTTAGCTTTGCATT
                            CT
                                          CTG
#>
              1
                             8
                                           3
#>
             CTGA
                          CTGAT
                                      CTGATTTT
#>
                                           5
              11
                            16
#>
         CTGATTTTT
                       CTGATTTTTA
                                    CTGATTTTTAG
#> CTGATTTTTAGATTTGCATT
                     CTGATTTTTAGC
                                   CTGATTTTTAGCT
#>
#>
      CTGATTTTTAGCTT
                   CTGATTTTTAGCTTT
                                 CTGATTTTTAGCTTTG
#>
              10
                                           4
                 CTGATTTTTAGCTTTGCA CTGATTTTTAGCTTTGCAT
#>
    CTGATTTTTAGCTTTGC
#> CTGATTTTTAGCTTTGCATT CTGATTTTTAGCTTTGCAAT
                                 CTGATTTTTAGCTTTT
                                           1
#> CTGATTTTTAGCTTTTCATA
                      CTGATTTTTAT
                            1
#> CTGATTTTTATCTTTGCATT CTGATTTTTATCTTTTGATT CTGTGTTTTAGATTTGCATT
                             1
                                           1
                          CTGAA CTTATTTTTATCTTTGCATT
#> CTGTTTTTTTTTTTTGCAAT
              1
#>
         TTAGCTTTG
H3122_results$breakpoint
#> break_pos
#> 42299750 42299757
```

```
#> 202 1

H3122_results$read_depth

#> [1] 251

HCC827_results

#> [1] "No EML4-ALK was detected"
```